

The universal pro - the ACROBAT 3000 from Ondal.

One system for all applications - the ACROBAT 3000 from Ondal is a reliable, universal system that provides for a broad range of options for medical applications. This wall-mounted or ceiling-mounted system supports monitors, camera systems and surgical lights in clinical applications.



ACROBAT 3000:

- Designed for loads up to 30 kg
- more than 360-degree swivel angle of support arms
- 330-degree swivel angle of support arms with anti-twist safeguard due to interior supply lines
- three-, five-, seven- or nine-pole commutators in swivel joints for power supply and signal transmission

Advantages of ACROBAT 3000:

- High-quality bearings for ease of movement.
- The precision parallel guide during height adjustment ensures the optimal positioning of terminal equipment.
- Interior cabling reduces potential trip hazards and lets you clean the system more thoroughly.

ondal[®]
Medical Systems

Movement – light as a feather