


## 31.5" 4K model FS-P3101D

FSN Medical Technologies has solutions for managing your operating room video signals. Our products are engineered and built for compatibility with other highly specialized surgical and diagnostic equipment used in surgical suites, operating rooms, emergency rooms, and procedural facilities. FS-P3101D, our 31" 4K medical display, is uniquely equipped to handle tasks in the demanding OR environment.



- 4K resolution (4096 x 2160 pixels).
- High brightness LED technology.
- DVI and DisplayPort inputs.
- Rapid signal detection, robust mode tables.
- Built to last, rugged aluminum housing.
- Calibrated to clinical color.
- Image mirror, rotate, pan, zoom, freeze, picture-in-picture.



### General Features

Panel	31 inch TFT LCD (LED)
Resolution	4096 x 2160 pixel
Power Supply	AC/DC Adaptor (AC 100-240V~, DC 24V 6.25A)
Control Key	Power, Menu, PIP, Up, Down, Plus, Minus, Input
Unit Dimension	761.5(W) x 442.0(H) x 68.0(D) mm, 29.98(W) x 17.40(H) x 2.67(D) inch
Average Brightness	350 cd/m2
Response Time (G-to-G)	20 ms
Viewing Angle (CR>10)	R/L 178, U/D 178
Weight	12.2 kg, 26.89 lbs.
Contrast Ratio	1500 : 1 typical
Pixel Pitch	0.1704 x 0.1704

