


MONITEUR 31.5" 4K modèle FS-P3101D

FSN propose des solutions pour gérer vos signaux vidéo dans la salle d'opération. Les produits FSN sont conçus et construits afin d'être compatibles avec d'autres équipements de haute technologie en chirurgie et diagnostics utilisés au bloc opératoire, dans les salles d'urgence et dans les salles d'endoscopie. Le moniteur FS-P3101D, avec ses 31,5" et sa résolution d'image médicale 4K, est particulièrement bien équipé pour gérer les tâches dans un environnement exigeant :

- Résolution 4K (4096 x 2160 pixels).
- Technologie LED haute luminosité.
- Entrées DVI et DisplayPort
- Détection de signal rapide
- Construit en aluminium pour durer
- Calibrage de la couleur pour l'environnement clinique
- Image miroir, rotation, vision panoramique, zoom, gel sur image, image dans l'image.



General Features

Panel	31 inch TFT LCD (LED)
Resolution	4096 x 2160 pixel
Power Supply	AC/DC Adaptor (AC 100-240V~, DC 24V 6.25A)
Control Key	Power, Menu, PIP, Up, Down, Plus, Minus, Input
Unit Dimension	761.5(W) x 442.0(H) x 68.0(D) mm, 29.98(W) x 17.40(H) x 2.67(D) inch
Average Brightness	350 cd/m2
Response Time (G-to-G)	20 ms
Viewing Angle (CR>10)	R/L 178, U/D 178
Weight	12.2 kg, 26.89 lbs.
Contrast Ratio	1500 : 1 typical
Pixel Pitch	0.1704 x 0.1704

