

Monitor Support

YK500

Medical  PM™

The YK500 has been designed to allow video screen upgrades; it considerably reduces the migration cost should you purchase new video equipment or upgrade to new monitors or articulated arms. YK500 minimizes installation time which eliminates the need of OR and Endoscopy shut down.

Specifically created to assist monitor upgrades, the YK500 spring arm interface is adaptable to most available spring arms on the market. It also avoids the replacement of the entire monitor support when comes time to install a new monitor flat panel.



YK500 is compatible with the following articulated arms:

- Ondal™
- Stryker™
- Dr Mach™
- Berchtold™
- Dräger™
- Heraeus™
- ALM Medical™

YK500 for which monitors ?

Compatible with 23", 26" and 32" monitors, the **YK500** can adapt to many screen dimensions thanks to its unique hanging feature, fitting on several models.



Monitor Support

YK500

Features

Hidden cabling in the **YK500** arm, preventing accidental shutdowns during procedures, maintenance or cleaning

Sterile handle adaptor

Sterile handle (not included)

Back splash that complies with infection control regulations and that avoids damaging the cables when cleaning (tailored covers also available)

Adjustable stop according to users' needs in 90 degree intervals

Assisted tilting motion

It's all in the handle...

The **YK500** is equipped with an adaptor for sterile handle that can be installed to the right or to the left of the monitor, according to user's preference.

Thanks to our adaptor, this handle can fit the same sterile handle that is used by most light manufacturers. This asset avoids the use of different sterile handles in the rooms, eases the sterilization process as well as simplifies stock inventory and manipulations.

Should you install a new surgical light, the **YK500** part receiving the sterile handle can be replaced to fit with the new light sterile handle.

We are developing a sterile handle adaptor fitting with the following surgical lights:

- Stryker™
- Berchtold™
- Heraeus™
- Steris™
- Dr Mach™
- ALM Medical™



Included with the YK500:

- Spring arm interface
- Yoke YK500
- Sterile handle adaptor extension
- Heraeus™ sterile handle adaptor
- Splash guard

Optional:

- Sterile handle adaptor for other brands than Heraeus™ surgical lights
- Tailored cover according to flat panels



Medical PM™

Toll Free: 1 (855) 404-7676

Email: info@medical-pm.com
www.medical-pm.com